Hera

Heraenium[®] CE, Heraenium[®] EH, Heraenium[®] NF, Heraenium[®] Laser

(DE) Gebrauchsanweisung	02	(CZ) Návod k použití	74
(GB) Instructions for use	08	(HU) Használati utasítás	80
(FR) Mode d'emploi	14	(LV) Lietošanas instrukcija	86
(ES) Instrucciones de uso	20	(LT) Naudojimo instrukcija	92
(IT) Istruzioni per l'uso	26	(PL) Instrukcja obsługi	98
(PT) Instruções de uso	32	(RU) Инструкция по применению	104
(NL) Gebruiksaanwijzing	38	(UA) Інструкція з використання	110
(SE) Bruksanvisning	44	(TR) Kullanma talimatı	116
(DK) Brugervejledning	50		
(NO) Bruksanvisning	56		
(FI) Käyttöohjeet	62		
(GR) Οδηγίες Χρήσης	68		

KULZER
MITSUBI CHEMICALS GROUP

Heraenium® CE, Heraenium® EH, Heraenium® NF, Heraenium® Laser

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Anwendung:

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Type 5 alloys (as per EN ISO 22674)

Usage:

Heraenium CE and EH/NF/Laser are cobalt chrome alloys for fabricating removable dentures. They exhibit high strength and outstanding elastic recovery. Heraenium CE and EH/NF/Laser are also highly elastic as well as being easily trimmed and welded.

1. Waxing up

CoCr frameworks are waxed up on an investment model according to the usual anatomical and mechanical criteria for CoCr dentures. Special KF casting wax is recommended for waxing and sprueing. KF casting wax is infiltrated with a grain refiner that greatly improves the structure of the casting and also the mechanical properties of the cast alloys.

2. Sprueing

Only two sprues are required with any type of framework. The diameter of the sprues should be 3.5 mm and the length 20–35 mm.

Feeder sprues with reservoirs are added to the thicker sections – the reservoirs should be 5.4 or 5.9 mm in diameter and the sprues 2.5 or 2.9 mm in diameter. Ventilation sprues, 0.8 mm in diameter and 1 cm in length, lead to a pressure relief sprue with a diameter of 1.2 mm, which runs around the framework. It is essential that the transitions from the sprue former to the sprues and from the sprues to the pattern are smooth. When attaching the sprues, ensure that the tip of the sprue former is approx. 3–4 mm above the highest point of the wax pattern.

3. Investing and preheating

The mixing ratios, expansion control and preheating schedules etc. vary depending on the investment material. The data for the investment materials are shown in the relevant instructions.

4. Melting and casting with vacuum induction casting machines

The Kulzer induction heated, vacuum pressure casting machines with ceramic crucibles for non-precious metal alloys are most suitable for melting and casting CoCr alloys.

	Heraenium CE	Heraenium EH	Heraenium Laser	Heraenium NF
1. Pre-melting	Pre-melting is discontinued when the last ingot is half melted			Melt metal fully (no visible contours)
2. Start casting	Interrupt pre-melting Open the casting drum Insert the mould Commence main melt			
3. Start casting by hand	Once the shadow disappears	3 seconds after the shadow disappears		6 seconds after the shadow disappears*

Important! Set the vacuum to 250 mbar. The vacuum must not be set higher!

* If the oxide crust splits during pre-melting or main melting, this need not influence your evaluation of the melting process

5. Trimming the casting

Unfinished castings are sandblasted using special abrasive M or 250 μm Al_2O_3 abrasive at a maximum pressure of 4.0 bar. This ensures that the optimum edge definition and surface structure of the casting remains intact. Metal dust is hazardous to the health. For finishing and sandblasting use a suitable exhaust system and/or a type FFP3-EN 149-2001 dust mask.

6. Soldering

When soldering Heraenium cast units to Heraenium or wire clasps, the units can first be secured in position by spot welding and then soldered using cobalt chrome/gold solder 910 and Kulzer Hera SLP 99 flux.

Heraenium cast units are soldered to precious dental alloy units with a flame using cobalt chrome/gold solder 750 and Hera SLP 99 flux. Pre-soldering is unnecessary. When soldering high-palladium-content bonding alloys to Heraenium, the solder junction of the bonding alloy must be pre-soldered with Herador 1060 before porcelain firing and post-soldered with cobalt chrome/gold solder 750.

7. Laser welding

The power of the pulse should be set according to the thickness of the units. Prolonging the pulse duration increases the volume of metal melted. This prolongs the solidification times, which has a positive effect in the welding seam. Before laser welding, all traces of solder must be completely removed from the seam area. Palladium-based alloys can never be welded to CoCr dentures.

Recommended settings (Herapuls):	Focal diameter:	\varnothing 1,1–1,3 mm
	Pulse duration:	12–15 ms
	Power duration:	depends on the thickness of the material

7a. Laser welding Heraenium CE, Heraenium EH, Heraenium NF

Please note! The carbon in the alloys may cause formation of carbide in the welding seam. This may cause the welding seam to become brittle. To prevent this, always weld with Kulzer special laser welding rod.

7b. Laser welding with Heraenium Laser

Please note! Welding without additives is possible by direct contact with joint areas. If there is a gap, always weld with Kulzer special laser welding rod or with similar added material.

8. Using casting buttons

Heraenium alloys must only be cast once. Kulzer repurchases casting buttons and sprues. Please ask for the current prices.

9. Contraindications, side effects and interactions with other dental alloys

In cases of hypersensitivity (allergy) to the constituents of Heraenium alloys, do not continue their use.

In individual cases, hypersensitivity reactions (allergies) and electrochemically-induced local dysaesthesia have been reported such as changes in taste and irritation of the oral mucosa.

Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Chemical composition								
	Chemical composition in % by weight							
	Co	Cr	Mo	Mn	Si	C	N	Ta
Heraenium CE	63,5	27,8	6,6	0,6	1,0	0,3	0,2	–
Heraenium EH	63,5	28,0	6,5	0,6	1,0	0,15	0,25	–
Heraenium NF	63,35	29,0	5,0	0,6	1,0	0,25	0,3	0,5
Heraenium Laser	63,5	28,0	6,5	0,6	1,05	< 0,05	0,3	–

Technical data							
	Density	Melting range °C	Casting temp. °C	Hardness	0,2 % yield strength	Elongation at rupture %	Young's modulus GPa
	g/cm ³			HV 10			
Heraenium CE	8,0	1330-1380	1530	380	580	4	228
Heraenium EH	8,0	1330-1380	1530	310	620	7	228
Heraenium NF	8,1	1330-1380	1530	360	720	5	230
Heraenium Laser	8,0	1330-1380	1530	340	610	12	220

Preheating temperature for all alloys 950 °C – 1030 °C **

** Depending on casting system.

10. Notes

The CoCr-based alloys may contain traces of nickel (<0.1 % be mass).
Please consult the material safety data sheets or national regulations for disposal.

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Heraenium® = registered trademark of Kulzer GmbH
Heraenium® = Kulzer GmbH kuruluşunun tescilli ticari markasıdır.
Kulzer GmbH