pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque

REF 032 Rev 011 Print Date: 30.10.2017

INTRODUCCIÓN
 All restorations made of pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque must go through a final sintering cycle prior to finishing. The sinter firing must only be conducted in high-temperature ovens authorized for this purpose. You may use the firing process specified by pri®dentita®. Please also note the oven manufacturer’s specifications. All standard programmable dental laboratory ovens may be used.

Sintering ovens vary in their performance, so it is highly recommended we use calibrate periodically to ensure the recommended cycle is carried out correctly. Follow the manufacturer’s recommended calibration instructions.

WARNING: Sintering furnaces must be located in a fireproof well-ventilated area. Do not open the oven or remove the sintered restoration before the furnace has cooled down sufficiently for safe handling to avoid the risk of burns.

Recommended firing program:
 From room temperature to 1.450 °C at 10 °C/min, hold time of 2 hours at firing temperature, allow to cool to room temperature at 10 °C/min. Natural cooling begins from approx. 600 °C depending on the oven used.

Rapid sintering (optional):
 Rapid sintering may be conducted subject to the following conditions:
 • cross-sectioned up to 3 units from room temperature to 1.500 °C at 10 °C/min, hold time of 30 minutes, allow to cool to room temperature at 40 °C/min.

WARNING: Always comply with the above sinter cycles, otherwise the material may become weakened and may break in the patient’s oral cavity.

You are strongly advised not to use coloring liquids as this can have a negative effect on translucency or shading.

If restorations have been fabricated in the sintering oven using liquid stains, a cleaning firing with diamond tipped grinding powder is recommended before pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque is sintered in the same oven.

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This device is intended to be sold only for processing by trained persons, milling laboratories or centers authorized to process dental prostheses.

INTRODUCTION
 Please read these instructions completely and thoroughly before using the product and observe them at all times. Improper use of the product and failure to observe the information provided may affect the quality and reduce the lifetime of the dental prosthesis.

The product must be used in accordance with the instructions and only for the specified purpose for which it was designed. pri®dentita® GmbH shall not be liable for any consequential damages or injury arising from any use or misuse of this product. In using pri®dentita® GmbH’s product you assume the responsibility of ownership and operation, and you agree to hold pri®dentita® GmbH harmless for any incident or action involving the use of pri®dentita® GmbH’s product.

Please keep these instructions in a safe location for further reference over the entire lifetime of the product. Pass on this information to any future owner, processor or user of this product or of other processed products made from it.

Please be aware of the different types of risk involved in the use of the product.
 • **WARNING:** Indicated a hazardous situation that, if not avoided, may result in serious injury.
 • **CAUTION:** Indicates a hazardous situation that, if not avoided, could result in minor or moderate injury or in property damage.

PRODUCT DETAILS
 Round blank, diameter: 98,5 mm, available with or without step
 Shades: white, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4
 Heights: 10, 12, 14, 16, 18, 20, 25 mm
 multicolor: A light (shade range: A1-A3), B light (shade range: B1-B2), C light (shade range: C1-C2), D light (shade range: D2-D4), A dark (shade range: A3.5-A4), B dark (shade range: B3-B4), C dark (shade range: C3-C4)
 Heights: 14, 18, 20, 25 mm

CERAMIC VENEERS
 All veneer ceramics recommended for zirconium oxide ceramics may be used.

FIXING
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque restorations may be cemented conventionally with a zinc phosphate or glass ionomer cement or with self-adhesive fixing and full coverage material. Ensure that surfaces are retentive and that stump height is at least 3 mm.

WARNING: Do not use a final restoration if there are any cracked, fractured or broken areas or discolorations visible on the product. If any defect is identified, the framework must not be used on a patient. Otherwise you will risk injuries in the patient’s oral cavity or the inhalation of the product or parts of it.

CONTRAINDICATIONS
 • Do not use the restoration
 • for the fabrication of implants
 • in patients with parafunctional habits (e.g. bruxism)
 • in patients with a known intolerance to components
 • with inadequate preparation
 • with insufficient space available in the patient’s mouth
 • in patients with inadequate oral hygiene

Otherwise you will risk severe patients’ injuries.

TECHNICAL SUPPORT
 For technical support contact manufacturer.

Manufacturer: pri®dentita® GmbH Meisenweg 37 70771 Leinfelden-Echterdingen Germany

Phone +49 (0) 711 320 656 0 Fax +49 (0) 711 320 656 99 E-Mail info@pri®dentita.com

REF 032 Rev 011 Print Date: 30.10.2017

INDICACIONES FOR USE
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque are pre-sintered zirconium oxide blanks for use in CNC (Computer Numerical Controls) milling machines to fabricate crowns, bridges, inlays and onlays for dental prostheses for restoration purposes only.

GENERAL INFORMATION
 Check the delivery immediately on receipt to ensure that:
 • all components are present
 • the packaging and the product are intact

WARNING: Do not use if there are any cracked, fractured or broken areas or discolorations visible on the product. If any defect is identified, the blank must not be used for fabrication of dental prostheses. Processing a cracked, fractured, broken or discolored product may result in a defective restoration with a risk of injury for patients.

If you detect any defect in the product, please contact your dealer or the manufacturer. Should there arise particular problems which are not covered sufficiently by these instructions, the matter should be referred to the manufacturer.

STORAGE CONDITIONS
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque blanks should be stored in their original packaging. Ensure that:
 • they are stored in a dry place
 • temperatures are between 5 °C (41 °F) and 50 °C (122 °F)
 • they are not subject to excessive vibration

CAUTION: Do not store in a damp area. Damppness may affect the product. Do not store near petrol contamination sources as they may affect the product.

MATERIAL PROPERTIES
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque round blanks meet the following specifications at completion of the high-density sintering process

Material characteristics:
 Flexural strength: > 1150 MPa 10⁻³ K⁻¹
 CTE: 10⁻⁶ K⁻¹
 Chemical composition: Zirconium dioxide (ZrO₂/HfO₂): 94,1 – 94,65 % Yttrium oxide (Y₂O₃): 4,65 – 5,95 % Aluminium oxide (Al₂O₃): 0,4 % Oxide other: < 0,7 %

The proportions of the single components in the total quantity of the components in a blank may vary within the ranges described above; however, the total quantity of components in each single blank will not exceed 100%.

MACHINING
 As pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque is made of a sensitive high performance material always handle it with care. Avoid handling with wet hands. This product must only be managed by trained technicians. The safety warnings in these instructions for use must be complied with to process bare sole responsibility for the use of the product. The manufacturer accepts no liability for defective results since it has no influence on the machining process.

CONSTRUCTION
 Systems must be constructed of the following specifications:

	Crowns	Mayland-bridges	3-unit bridges
Minimum framework thickness	anterior 0,4 mm posterior 0,6 mm	0,4 mm 0,6 mm	0,6 mm 0,6 mm
Connectors	anterior – 6 mm ² posterior 9 mm ²	6 mm ² 9 mm ²	6 mm ² 9 mm ²
Framework design	Anatomically shaped (supporting veneer ceramics); fully anatomical		

For technical support contact manufacturer.

Manufacturer: pri®dentita® GmbH Meisenweg 37 70771 Leinfelden-Echterdingen Germany

Phone +49 (0) 711 320 656 0 Fax +49 (0) 711 320 656 99 E-Mail info@pri®dentita.com

REF 032 Rev 011 Print Date: 30.10.2017

INDICACIONES FOR USE
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque are pre-sintered zirconium oxide blanks for use in CNC (Computer Numerical Controls) milling machines to fabricate crowns, bridges, inlays and onlays for dental prostheses for restoration purposes only.

GENERAL INFORMATION
 Check the delivery immediately on receipt to ensure that:
 • all components are present
 • the packaging and the product are intact

WARNING: Do not use if there are any cracked, fractured or broken areas or discolorations visible on the product. If any defect is identified, the blank must not be used for fabrication of dental prostheses. Processing a cracked, fractured, broken or discolored product may result in a defective restoration with a risk of injury for patients.

If you detect any defect in the product, please contact your dealer or the manufacturer. Should there arise particular problems which are not covered sufficiently by these instructions, the matter should be referred to the manufacturer.

STORAGE CONDITIONS
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque blanks should be stored in their original packaging. Ensure that:
 • they are stored in a dry place
 • temperatures are between 5 °C (41 °F) and 50 °C (122 °F)
 • they are not subject to excessive vibration

CAUTION: Do not store in a damp area. Damppness may affect the product. Do not store near petrol contamination sources as they may affect the product.

MATERIAL PROPERTIES
 pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque round blanks meet the following specifications at completion of the high-density sintering process

Material characteristics:
 Flexural strength: > 1150 MPa 10⁻³ K⁻¹
 CTE: 10⁻⁶ K⁻¹
 Chemical composition: Zirconium dioxide (ZrO₂/HfO₂): 94,1 – 94,65 % Yttrium oxide (Y₂O₃): 4,65 – 5,95 % Aluminium oxide (Al₂O₃): 0,4 % Oxide other: < 0,7 %

The proportions of the single components in the total quantity of the components in a blank may vary within the ranges described above; however, the total quantity of components in each single blank will not exceed 100%.

MACHINING
 As pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque is made of a sensitive high performance material always handle it with care. Avoid handling with wet hands. This product must only be managed by trained technicians. The safety warnings in these instructions for use must be complied with to process bare sole responsibility for the use of the product. The manufacturer accepts no liability for defective results since it has no influence on the machining process.

CONSTRUCTION
 Systems must be constructed of the following specifications:

	Crowns	Mayland-bridges	3-unit bridges
Minimum framework thickness	anterior 0,4 mm posterior 0,6 mm	0,4 mm 0,6 mm	0,6 mm 0,6 mm
Connectors	anterior – 6 mm ² posterior 9 mm ²	6 mm ² 9 mm ²	6 mm ² 9 mm ²
Framework design	Anatomically shaped (supporting veneer ceramics); fully anatomical		

For technical support contact manufacturer.

Manufacturer: pri®dentita® GmbH Meisenweg 37 70771 Leinfelden-Echterdingen Germany

Phone +49 (0) 711 320 656 0 Fax +49 (0) 711 320 656 99 E-Mail info@pri®dentita.com

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ADJUSTS
 If it was necessary to make adjustments in the restorations densely sintered, these adjustments can only be made using machine tools with the aim of avoiding damage to the materials by localized overheating or excessive force on the surface of the restoration. No use any filing tools, as they will damage the surface of the restoration.

Basic rules for handling sintered materials:
 • work at low pressure.
 • use diamond-tipped grinding tools in good condition only.
 • avoid sharp edges at the finished product to prevent patient injuries.
 • do not machine interdental connectors.
 • prior to cementing the restoration in place, clean it and sandblast the internal surfaces in accordance with the following parameters: jet pressure 1 bar, particle size < 50 µm and distance approx. 10 mm.

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FLUENING
 The restorations pri® multidis ZrO₂ Translucet / pri® multidis ZrO₂ Opaque can be sintered in a furnace with a zirconium oxide axis, but only with a self-adhesive and fully adhesive luting agent. Ensure that the surfaces have a good fixation and that the stump height is at least 3 mm.

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