Signum® composite / composite flow / matrix®

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Signum® composite / composite flow / matrix®

The products of the brand Signum are light curing veneering composites for crowns and bridges.

Fields of Application

- Full and partial veneering for fixed metal-supported crowns and bridges.
- Veneering of resin based frames (PMMA).
- Veneering of removable prostheses (telescope/coronic crowns and attachments).
- For the veneering of Implant supra-/tertiary structures.
- Zirconia-supported longterm provisionals.
- Temporary anterior, posterior crowns (w/o frame).
- Characterisation of artificial teeth.

Note

The metal supported veneering is permitted up to the second premolar teeth without limitation, unless the molars are expected to be subject to undue stress. Signum composite, Signum composite flow and Signum matrix may also be applied for the occlusal veneering of molars to obtain aesthetically pleasing results. In case of any foreseeable physiological strain or forces acting on the masticatory surfaces of the molars (e.g. bruxism or parafunctions), veneer contact must be ruled out or metal stops applied. Fully veneered pontics should be separated interdentially prior or after the polymerisation in order to minimise the risk of stress cracks. Especially for the press to – or injection technique a sufficient interdental separation for long-span bridges must exist!

Composition

Signum composite

Signum dentine, Signum margin, Signum enamel, Signum value, Signum effect: Multifunctional methacrylic acid esters, photo initiators, stabilising agents, inorganic pigments (cadmium free).

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Fillers: 74% w/w silicon dioxide, silicon dioxide rheologically active type, splitter pre-polymerised.

Signum composite flow / Signum matrix

Signum mameion dentine, Signum secondary dentine, Signum opal incisal, Signum opal transparent, Signum composite flow: Multifunctional methacrylic acid esters, photo initiators, stabilising agents, inorganic pigments (cadmium free).

Fillers: 64% w/w, silicon dioxide, Ba-Al-B-silicate glass.

Framework Conditioning

In conjunction with 0.4 mm retention beads (please refer to detailed information at the Signum instructions for use) use the veneering composites of Signum only together with approved bonding systems, (Signum® metal bond, Signum® zirconia bond). Follow the Instructions of use of the individual bonding systems.

Hazard/Safety information

Signum composite / matrix

May cause an allergic skin reaction. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. Wash contaminated clothing before reuse. If skin irritation or rash occurs: Get medical advice/attention.

Signum composite flow

Can cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Use personal protective equipment as required. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Framework Modelling, Material Thickness, Finishing and Polishing

Please refer to detailed information at the Signum instructions for use.

Storage

Do not store material above 25°C (77°F) in original containers. This also applies to containers already opened and resealed. The products of Signum must not be

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