Keramik-Rohlinge
CAM-Blanks
Grezzi in ceramica
Piezas brutas cerámicas
Blocs pour technologie FAO
CAM-Blanks
Керамические заготовки

Metoxit AG www.metoxit.com
CH-8240 Thayngen – Switzerland

high tech ceramics

Made in Switzerland

high tech ceramics

high tech ceramics

high tech ceramics

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GB Instruction for use Metoxit Z-CAD® CAM-Blanks

- Fix the blank according to the CAM-System instructions, set the enlargement-factor.
- Start the machining.
- After machining, dismount the blank with the framework.
- Remove the framework from the blank with a diamond disc.
- De-dust the framework with an air-gun.
- Dry the framework completely before firing.
- Position the framework in a sinter-bed inside the crucible.
- Programme your furnace according to its manual and Metoxit sinter instructions.
- Sinter the framework.
- Wait until cool down is completed before removing the sintered framework.
- Inspect the framework in respect of flaws, wall thickness, seat and margin fit. If necessary, make small adjustments with a water-cooled diamond-tool.
- Rinse the framework with water and dry it. It is now ready for finishing or veneering!

Safety measures

- Design the framework to maximize cross-sections.
- Carefully read the MSDS in the download area: www.metoxit.ch
- Carefully read the Declaration of conformity: http://www.metoxit.com/downloads/zertifikat/Konformitaetserklaerung_Z-A-CAD_Flow_Pen_Liquid_2017.pdf

Sinter instructions

		Ramp Rate °C / h	T _{max} °C	Dwell Time h	Cooling
Z-CAD® Standard	Crowns	600	1500	1	Natural furnace cooling
			1450	2	
	Bridges	300	1500	1	
			1450	2	
Z-CAD® HTL / FTL / smile high translucent	Crowns	600	1450	2	
	Bridges	300	1450	2	

For Speed-Sintering a furnace specific trial should be conducted.

Indications for use

Metoxit Z-CAD®-blanks are indicated for use as a substructure for porcelain fused ceramic fixed dental restorations. Limitations are listed in table 1. All white Z-CAD® Blanks are indicated to be shaded with Z-CAD® Liquid, see appropriate IFU for further information.

Table 1: Indications of use and maximum number of pontics.

Material	Single Crowns		Bridges, Number of Pontics	
	Anterior	Posterior	Anterior	Posterior
Z-CAD® Standard	X	X	2	2
Z-CAD® HTL / FTL	X	X	2	2
Z-CAD® smile	X	X	1 (max. 3 units)	not to be indicated

X indicated; digits show maximum number of pontics.

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Contraindications

Contraindications are all applications not indicated in table 1.

Z-CAD® and FlowPen® are registered names and trademarks of Metoxit AG.

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Z-CAD® y FlowPen® son nombres registrados y marcas comerciales de Metoxit AG.
Z-CAD® e FlowPen® são nomes e marcas registados da Metoxit AG.